



WE SHAPE THE METAL TO SERVE HUMANITY



ELECTROMED INSTRUMENTS (PVT) LTD.

MADE BY PROFESSIONAL'S FOR PROFESSIONAL'S



WELCOME TO **ELECTROMED INSTRUMENTS (PVT) LTD.**

MISSION STATEMENT:

Our Company **ELECTROMED INSTRUMENTS** is striving and will be to develop a selected range of “**ELECTROSURGERY**” and “**GENERAL SURGICAL**” instruments from time to time keeping in view the new development and procedures in the medical world. New developments in surgical and General Surgical instrumentation introduced that will distinguish “**ELECTROMED INSTRUMENTS**” as a one-step forward company in Pakistan.

Our management will concentrate on manufacturing instruments in-house with the latest technology. **ELECTROMED INSTRUMENTS** will be introduced as a world-class brand under the international distributor’s network.

Quality management system together with international certifications will be updated as and when necessary. Quality assurance at all manufacturing steps will be the key to claiming our product reliability to consider “**ELECTROMED INSTRUMENTS**” as one of the best suppliers.

Our workers satisfaction at all management levels will be monitored under internal management system to enhance the productivity and skill level.

ELECTROMED INSTRUMENTS will achieve annual growth in manufacturing volumes and exports by setting up realistic targets. We pray for the prosperity and excellence of our management, workers and business tasks.

QUALITY POLICY:

Our Company believes in using “**TECHNOLOGY**” to produce instruments with “**PRECISION**” for the “**SATISFACTION**” of Doctors and Patients. Our products will be cost-effective with efficient delivery as desired by the Customers. After-sale Service and intelligent stock management will be monitored for continuous satisfaction.

TIME SAVING OPEN SURGERY:

Electro Surgical Instruments help Doctors to save time in open surgery with the best performance results. The Doctors are supported in a way that they perform to the best of their skills by using the latest techniques with our instruments.

ONE STEP COAGULATION AND PROCEDURES:

Electro Surgical Instruments reduce the usage of a number of instruments in surgery procedures. **ELECTROMED INSTRUMENTS** has brought prosperity in the medical field that Doctors can’t resist using our instrumentation.

PRIORITIZING PATIENT SAFETY:

The Bipolar instruments of Our Company intelligently control the current flow and reduce lateral spread wear-resistant ceramic inlay in a way that electrical cross circuits never confuse. Patient safety is our priority.

SURGEON’S COMFORT AND SATISFACTION:

Our Company has designed patterns that are acceptable to any doctor in any part of the world. Our designs bring a better grip and feel to the doctors and enhance their skill level which ultimately leads to their satisfaction. We offer ready-to-use products.

COST EFFECTIVENESS AND RELIABILITY:

Our goal is to continuously offer cost-effective products together with enhanced reliability all the times. Our materials and workmanship are of world-class standards. International brands and Doctors prefer to depend on **ELECTROMED INSTRUMENTS** over any other brand.



EMI-120
 50x
 1x
 Biopolar
 550



EMI-130
 50x
 1x
 Biopolar
 550



EMI-140
 50x
 1x
 Biopolar
 550



EMI-150
 50x
 1x
 Biopolar
 550



EMI-160
 50x
 1x
 Biopolar
 550



EMI-170
 50x
 1x
 Biopolar
 550



EMI-171
 50x
 1x
 Biopolar
 550



EMI-200
 50x
 1x
 Biopolar
 550



EMI-210
 50x
 1x
 Biopolar
 550



EMI-240
 50x
 1x
 Biopolar
 550

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGTH (CM)	ITEM CODE
American Bipolar Forceps Coaptation	0.3, 0.5	6, 8, 9, 10	11 CM	EMI-120
American Bipolar Forceps Iris (Straight)	0.3, 0.5	6, 8, 9, 10	11 CM	EMI-130
American Bipolar Forceps Iris (Cruved)	0.3, 0.5	6, 8, 9, 10	11 CM	EMI-140
American Bipolar Forceps Coaptation Nadler	0.3, 0.5	6, 8, 9, 10	12.5 CM	EMI-150
American Bipolar Forceps Mcpherson (Cruved)	0.3, 0.5	6, 8, 9, 10	10.5 CM	EMI-160
American Bipolar Forceps Mcpherson (Straight)	0.3, 0.5	6, 8, 9, 10	10.5 CM	EMI-170
American Bipolar Forceps Jeweller #7	0.3, 0.5	6, 8, 9, 10	13 CM	EMI-171
American Bipolar Forceps Jeweller #5	0.3, 0.5	6, 8, 9, 10	11.5 CM	EMI-200
American Bipolar Forceps Jewellers #4	0.3, 0.5	6, 8, 9, 10	11.5 CM	EMI-210
American Bipolar Forceps Adson	0.3, 0.5, 1.0, 1.5	6, 8, 9, 10	12.5 CM	EMI-240



Product Features:

Smooth surface, conductive and accurate coagulation can reduce the generation of eschar.

More Safety for high quality of insulated layer.

Different shape and size of forceps are suitable for various surgery.

Available in reusable and disposable type. Reusable forceps can withstand high pressure and steam sterilization.

Steam: Steam autoclave at 250°F / 121°C for 20 Minutes.

Flash: Steam autoclave at 275°F / 134°C for 10 Minutes.

Gas (ETO): Follow the manufacturer's guide for the unit which you are using.



EMI-250

50x
1x
Biopolar
550



EMI-260

50x
1x
Biopolar
550



EMI-270

50x
1x
Biopolar
550



EMI-280

50x
1x
Biopolar
550



EMI-290

50x
1x
Biopolar
550



EMI-300

50x
1x
Biopolar
550



EMI-340

50x
1x
Biopolar
550



EMI-360

50x
1x
Biopolar
550

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGTH (CM)	ITEM CODE
American Bipolar Forceps Semkins	0.5, 1.0, 1.5	6, 8, 9, 10	14 CM	EMI-250
American Bipolar Forceps Cushing (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-260
American Bipolar Forceps Gerld	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-270
American Bipolar Forceps Gerld (Curved)	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-280
American Bipolar Forcep Gerald Bayonet	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-290
American Bipolar Forceps YasarGil Bayonet	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-300
American Bipolar Forceps Jansen Bayonet	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-340
American Bipolar Forcep Scoville Greenwood	0.5, 1.0, 1.5	6, 8, 9, 10	19.5 CM	EMI-360



EMI-241-G.L.E.B.C
 50x
 1x
 Biopolar
 550



EMI-242-G.L.E.B.C
 50x
 1x
 Biopolar
 550



EMI-261-G.L.E.B.C
 50x
 1x
 Biopolar
 550



EMI-271-B.N.C
 50x
 1x
 Biopolar
 550



EMI-281-B.N.C
 50x
 1x
 Biopolar
 550



EMI-302-B.C
 50x
 1x
 Biopolar
 550



EMI-367
 50x
 1x
 Biopolar
 550



EMI-443-B.C
 50x
 1x
 Biopolar
 550



EMI-860-B.C
 50x
 1x
 Biopolar
 550



EMI-950-TW.B.C
 50x
 1x
 Biopolar
 550

DESCRIPTION	TIP WIDTH (MM)	TIPLength (MM)	ITEM LENGTH (CM)	ITEM CODE
American Bipolar Forceps (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	12 CM	EMI-241-G.L.E.B.C
American Bipolar Forceps (Curved)	0.5, 1.0, 1.5	6, 8, 9, 10	12 CM	EMI-242-G.L.E.B.C
American Bipolar Forceps (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-261-G.L.E.B.C
American Bipolar Forceps Gerld (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-271-B.N.C
American Bipolar Forceps Gerld (Curved)	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-281-B.N.C
American Bipolar Forceps Yasir Gill	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-302-B.C
TW American Bayonet Bipolar Forceps	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-367
American Adson Bipolar Forceps Non-Stick	0.5, 1.0, 1.5	6, 8, 9, 10	12.5 CM	EMI-443-B.C
European Bipolar Forceps (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-860-B.C
TW European Bayonet Forceps	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-950-TW.B.C



EMI-510
 50x
 1x
 Biopolar
 550



EMI-535
 50x
 1x
 Biopolar
 550



EMI-531
 50x
 1x
 Biopolar
 550



EMI-471
 50x
 1x
 Biopolar
 550



EMI-500
 50x
 1x
 Biopolar
 550



EMI-550
 50x
 1x
 Biopolar
 550



EMI-570
 50x
 1x
 Biopolar
 550



EMI-582
 50x
 1x
 Biopolar
 550



EMI-345
 50x
 1x
 Biopolar
 550



EMI-595
 50x
 1x
 Biopolar
 550

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGHT (CM)	ITEM CODE
American Bipolar Forceps Mcpherson (Straight)	0.3, 0.5, 1.0	6, 8, 9, 10	10.5 CM	EMI-510
American Bipolar Forcep Jeweller #5 Non Stick	0.3, 0.5, 1.0	6, 8, 9, 10	11.5 CM	EMI-535
American Bipolar Forcep Jeweller #5 Non Stick (Curved)	0.3, 0.5, 1.0	6, 8, 9, 10	11 CM	EMI-531
American Bipolar forceps Adson Two Pin Connector Non-Stick	0.3, 0.5, 1.0	6, 8, 9, 10	12 CM	EMI-471
American Bipolar Forcep Cushing Non Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	18 CM	EMI-500
American Bipolar Forceps Bayonet YasarGil Non-Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-550
American Bipolar Forcep Gerald Non Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	18.5 CM	EMI-570
American Bipolar Forceps Jenson Bayonet Non Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-582
American Bipolar Forceps New Bayonet	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	16.5 CM	EMI-345
A.C. American Boyonet Bipolar Forceps Three Holes Non-Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	24 CM	EMI-595



EMI-793
 50x
 1x
 Biopolar
 550



EMI-794
 50x
 1x
 Biopolar
 550



EMI-795
 50x
 1x
 Biopolar
 550



EMI-796
 50x
 1x
 Biopolar
 550



EMI-797
 50x
 1x
 Biopolar
 550



EMI-510-S.C
 50x
 1x
 Biopolar
 550



EMI-531-S.C.O
 50x
 1x
 Biopolar
 550



EMI-472-S.C
 50x
 1x
 Biopolar
 550



EMI-503-S.C.O
 50x
 1x
 Biopolar
 550



EMI-1105-N.R
 50x
 1x
 Biopolar
 550

DESCRIPTION	TIPWIDTH (MM)	TIPLENGTH (MM)	ITEM LENGTH (CM)	ITEM CODE
European L.C Jeweller #5 Bipolar Forceps	0.3, 0.5, 1.0	6, 8, 9, 10	12 CM	EMI-793
European L.C Jeweller #5 Bipolar Forceps (Curved)	0.3, 0.5, 1.0	6, 8, 9, 10	12 CM	EMI-794
European L.C with 3 Holes Bipolar Forcep (Straight)	0.3, 0.5, 1.0	6, 8, 9, 10	16 CM	EMI-795
European L.C Bayonet with 3 Holes Bipolar Forcep	0.5, 1.0, 1.5	6, 8, 9, 10	16 CM	EMI-796
European L.C Bayonet with 3 Holes Bipolar Forcep	0.5, 1.0, 1.5	6, 8, 9, 10	19.5 CM	EMI-797
American Bipolar Forceps Mcpherson Non-Stick (Straight)	0.3, 0.5, 1.0	6, 8, 9, 10	10.5 CM	EMI-510-S.C
American Bipolar Forcep Jeweller #5 Non-Stick (Straight)	0.3, 0.5, 1.0	6, 8, 9, 10	11.5 CM	EMI-531-S.C.O
American Adson Bipolar Forceps with Five Holes Non-Stick	0.5, 1.0, 1.5	6, 8, 9, 10	12 CM	EMI-472-S.C
European Bipolar Forcep Non-Stick (Straight)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-1105-N.R



EMI-790

50x
1x
Biopolar
550



EMI-740

50x
1x
Biopolar
550



EMI-750

50x
1x
Biopolar
550



EMI-770

50x
1x
Biopolar
550



EMI-780

50x
1x
Biopolar
550



EMI-840

50x
1x
Biopolar
550



EMI-830

50x
1x
Biopolar
550



EMI-872-N

50x
1x
Biopolar
550



EMI-975-S.S

50x
1x
Biopolar
550



EMI-976

50x
1x
Biopolar
550

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGTH (CM)	ITEM CODE
European Bipolar Forcep Jeweller #5	0.3, 0.5, 1.0	6, 8, 9, 10	11.5 CM	EMI-790
European Bipolar Forceps Jewellery #5 (Curved)	0.3, 0.5, 1.0	6, 8, 9, 10	11.5 CM	EMI-740
European Bipolar Forcep Jeweller #7	0.3, 0.5, 1.0	6, 8, 9, 10	11 CM	EMI-750
European Bipolar Forcep Mcpherson	0.3, 0.5, 1.0	6, 8, 9, 10	10 CM	EMI-770
European Bipolar Forcep Adson	0.3, 0.5, 1.0	6, 8, 9, 10	12 CM	EMI-780
European Bipolar Forceps (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	15 CM	EMI-840
European Bipolar Forceps (Curved)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-830
European Bipolar Forcep (Straight)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-872-N
European Bipolar Forceps Bayonet	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-975-S.S
European Bipolar Forceps Bayonet (Curved)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-976



EUROPEAN NON-STICK BIPOLAR FORCEPS

NON-STICK Forceps

The newest NON-STICK Forceps extend the existing portfolio of NON-Stick forceps with a wide range of different lengths and shapes of tatteded and fine coagulation.

with proven NON-Stick technology based on using a special alloy in the tips, they guarantee excellent heat dissipation and reduce adhesion of tissue to a minimum.



EMI-1020

50x
1x
Biopolar
550



EMI-990

50x
1x
Biopolar
550



EMI-1105-N

50x
1x
Biopolar
550



EMI-1080

50x
1x
Biopolar
550



EMI-1101

50x
1x
Biopolar
550



EMI-1050

50x
1x
Biopolar
550



EMI-1051-N

50x
1x
Biopolar
550



EMI-1061-S.S

50x
1x
Biopolar
550



EMI-1062-S.S

50x
1x
Biopolar
550

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGTH (CM)	ITEM CODE
European Bipolar Forcep Jeweller # 5 Non Stick	0.3, 0.5, 1.0	6, 8, 9, 10	11.5 CM	EMI-1020
European Bipolar Forcep Adson Non Stick	0.3, 0.5, 1.0	6, 8, 9, 10	12 CM	EMI-990
European Bipolar Forcep Non-Stick (Straight)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-1105-N
European Bipolar Forcep Non-Stick (Straight)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	18 CM	EMI-1080
European Bipolar Forcep Non-Stick Straight (Curved)	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	21.5 CM	EMI-1101
European Bayonet Bipolar Forcep Non Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-1050
European Bipolar Forcep Bayonet Non Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-1051-N
European Bipolar Forceps Bayonet Non-Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-1061-S.S
European Bipolar Forceps Bayonet Angled Down Non-Stick	0.5, 1.0, 1.5, 2.0	6, 8, 9, 10	20 CM	EMI-1062-S.S

DOUBLE IRRIGATION BIPOLAR FORCEPS

EMI-693



EMI-694



EMI-692



Product Features:

1. Irrigation device lowers heat in the surgical site and avoids thermal tissue damage.
2. Continuous dripping moisturizes the tissue and minimizes the danger of dessication.
3. No adhesion and no eschar formation ensures clean and clear visualization in the surgical field.
4. Suitable for neurosurgery, plastic surgery, ENT, Vascular surgery and hepatobiliary surgery.
5. Also available in NON-STICK Tip.

Convert any regular Bipolar forceps into NON-STICK Irrigation bipolar forcep use "NS" extension for conversion available.

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGHT (CM)	ITEM CODE
Double Irrigation Boynet with 3M Cable, With Drip Pipe	0.5, 0.7, 1.0	7, 8	19 CM	EMI-693
Double Irrigation Boynet with 3M Cable, With Drip Pipe	0.5, 0.7, 1.0	7, 8	23.5 CM	EMI-694
Double Irrigation (Straight) with 3M Cable, With Drip Pipe	0.5, 0.7, 1.0	7, 8	20 CM	EMI-692



SINGLE IRRIGATION BIPOLAR FORCEPS



EMI-651



EMI-671




EMI-691-B.C



EMI-600




EMI-1110




EMI-1120

MONOPOLAR FORCEPS



EMI-370
50x
1x
Monopolar
550



EMI-445
50x
1x
Monopolar
550

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGHT (CM)	ITEM CODE
American Irrigation Bipolar Forceps Adson	0.5, 1.0, 1.5	6, 8, 9, 10	14 CM	EMI-651
American Caushing irrigation Bipolar Forceps	0.5, 1.0, 1.5	6, 8, 9, 10	16.5 CM	EMI-671
American Irrigation Bipolar Forceps Jenson Bayonet	0.5, 1.0, 1.5	6, 8, 9, 10	19 CM	EMI-691-B.C
American Irrigation Bipolar Forceps Jenson Bayonet Non-Stick	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-600
European Irrigation Bipolar Forcep (Straight)	0.5, 1.0, 1.5	6, 8, 9, 10	18 CM	EMI-1110
European irrigation Bipolar Forcep Bayonet	0.5, 1.0, 1.5	6, 8, 9, 10	20 CM	EMI-1120


MONOPOLAR FORCEPS

European Monopolar Female 4.0mm Connector	0.5, 1.0, 1.5	6, 8, 9, 10	21.5 CM	EMI-370
Jenson Boynet Angled Up Monopolar Forceps Screw Fitting	0.5, 1.0, 1.5	6, 8, 9, 10	21 CM	EMI-445



LAPAROSCOPIC ELECTRODS

ALSO AVAILABLE IN
TAFOLON COATED



WITHOUT SUCTION



REUSABLE

DISPOSABLE

EMI-3110

EMI-3111

L-Hook Electrode (33, 36, 39 cm)



EMI-3112

EMI-3113

J-Hook Electrode (33, 36, 39 cm)



EMI-3114

EMI-3115

Needle Electrode (33, 36, 39 cm)



EMI-3116

EMI-3117

Spatula Electrode (33, 36, 39 cm)



EMI-3118

EMI-3119

Ball Electrode (33, 36, 39 cm)



WITH SUCTION



REUSABLE

DISPOSABLE

EMI-3120

EMI-3121

L-Hook Electrode (33, 36, 39 cm)



EMI-3122

EMI-3123

J-Hook Electrode (33, 36, 39 cm)



EMI-3124

EMI-3125

Needle Electrode (33, 36, 39 cm)



EMI-3126

EMI-3127

Spatula Electrode (33, 36, 39 cm)



EMI-3128


EMI-3129

Ball Electrode (33, 36, 39 cm)





BIPOLAR ARTERY SEALER



EMI-4006



EMI-4007



EMI-4008



EMI-4009

DESCRIPTION	TIP WIDTH (MM)	TIP LENGTH (MM)	ITEM LENGTH (CM)	ITEM CODE
Reusable Non-Stick Tip Artery Sealer Black Color	4.0	40	22 CM	EMI-4006
Reusable Non-Stick Tip Artery Sealer Red Color	4.0	40	22 CM	EMI-4007
Reusable Non-Stick Tip Artery Sealer Purple Color	3.5	35	19 CM	EMI-4008
Reusable Polished Tip Artery Sealer Blue Color	3.5	35	19 CM	EMI-4009



All 2.4mm needle electrodes are also available with 4.0mm shaft. The article number is supplemented with "4" after the first two letters.



- EMI-03234** Arthroscopic Blade Electrode Connector 2.4mm Length 7 cm D:2.5mmx7mm Black Molded Non Stick Coated
- EMI-03236** Arthroscopic Blade Electrode Connector 2.4mm Length 7 cm D:2.5mmx17mm Non Stick Coated Black Insulated
- EMI-152-510** 4 mm Circle Loop Electrode 2.4 mm connecting shaft Black Molded
- EMI-152-520** 6 mm Circle Loop Electrode 2.4 mm Connecting shaft Black Molded
- EMI-152-530** 8 mm Circle Loop Electrode 2.4 mm connecting shaft Black Molded
- EMI-152-540** 6 mm Diamond Loop Electrode 2.4 mm connecting shaft Black Molded
- EMI-152-550** Fine Needle Electrode Connector 2.4 mm Length 7cm D: 0.25mm x 15mm Black Molded
- EMI-1767** 3 mm Ball Electrode Black Molded Connector Shaft 2.4 mm Length 7 cm
- EMI-03250** Schuco Knife Electrode 2.5mm x 18mm Connector Shaft 2.4mm Length = 70mm Black Molded
- EMI-3238** Needle Electrode Connector 2.4 mm Length 7 cm Teflon Coated, D: 1.3mm x 25mm
- EMI-3240** Arthroscopic Needle Electrodes 2.4mm connector Tungsten Tip D:0.8mm x 18mm Length 5.5cm
- EMI-03224** Arthroscopic Electrodes 2.4mm connector / Tungsten tip Needle Electrode Straight 6 cm D:0.8mm x 15mm
- EMI-3225** Arthroscopic Electrodes 2.4mm connector Needle Electrode Angled 7 cm
- EMI-03241** Needle Electrode 2.4mm connector Length 5.5 cm D: 0.7mm x 17mm Blue Molded
- EMI-03223** Arthroscopic Electrodes 2.4mm connector / Tungsten tip Needle Electrode Straight 11.5 cm D:0.8mm x 15mm
- EMI-3239** Arthroscopic Needle Electrodes 2.4mm connector Length 13.5 cm D: 0.9mm x 13 mm Blue Molded
- EMI-03243** Needle Electrode 2.4mm connector Angled Length 14 cm D: 0.8mm x 10mm Blue Molded
- EMI-1741** 3 mm Ball Electrode Length 13.5 cm Connector 2.4 mm
- EMI-1763** 5mm Ball Electrode molded Connector Shaft 2.4 mm Length 13.5 cm
- EMI-1761** 4mm Ball Electrode Non Stick Teflon Coated Connector Shaft 2.4 mm Length 13.5 cm
- EMI-1762** 5mm Ball Electrode Non Stick Teflon Coated Connector Shaft 2.4 mm Length 13.5 cm
- EMI-03251** Loop Electrode D:10mm X 12mm Length 16cm 2.4mm connector, Black insulated
- EMI-03242** Needle Electrode 2.4mm connector Black insulated Length 15 cm D: 0.8mm x 18 mm Tungsten tip Angled
- EMI-3220** Arthroscopic Blade Electrodes Dagger Tip Non Stick Coated Length 18 cm D: 2.3mmx18mm



- EMI-1783** Conisation Electrode 25mm width (B) 20mm Depth (A) , connector 2.4mm working Length =95mm Total Length 14 cm Blue Molded and Blue coated
- EMI-1782** Conisation Electrode 25mm width (B) 15mm Depth (A) , connector 2.4mm working Length =95mm Total Length 14 cm Blue Molded and Blue coated
- EMI-1781** Conisation Electrode 25mm width (B) 10mm Depth (A) , connector 2.4mm working Length =95mm Total Length 14 cm Blue Molded and Blue coated
- EMI-03254** Round Loop electrode 5mm , connector 4mm Length 7 cm Blue Molded
- EMI-03255** Round Loop electrode 10mm , connector 4mm Length 7 cm Blue Molded
- EMI-4-1769** 4mm Ball Electrode Blue Molded Connector Shaft 4 mm Length 5.5 cm
- EMI-1769** Fine Needle Electrode Connector 2.4 mm Length 7cm D: 0.25mm x 15mm Black Molded
- EMI-1770** 6mm Ball Electrode Blue Molded Connector Shaft 4 mm Length 5.5 cm
- EMI-03239** Arthroscopic Blade Electrode Connector 4mm Length 5.5 cm D:3.3mmx25mm Blue Molded
- EMI-03247** knife Electrode connector 4 mm Length 5.5 cm D: 3.3mm x 25mm Blue Molded
- EMI-03240** Arthroscopic Needle Electrodes 2.4mm connector Tungsten Tip D:0.8mm x 18mm Length 5.5cm
- EMI-03245** knife Electrode connector 4 mm Length 14 cm D: 2.5mm x 12mm Blue Molded
- EMI-03252** V-Shape Loop electrode 10mm x 12mm , connector 4mm Length 170mm Blue Molded and blue insulated
- EMI-03253** V-Shape Loop electrode 15mm x 15mm , connector 4mm Length 170mm Blue Molded and blue insulated
- EMI-1765** 8mm Ball Electrode Molded Connector Shaft 4 mm Length 14 cm
- EMI-3235** Needle Electrode Connector 4mm Length 15 cm D: 1.3mm x 20 mm Blue Molded
- EMI-03244** knife Electrode Angled connector 4mm Length 16 cm D: 2.5mm x 12mm Blue Molded
- EMI-03233** 5mm Ball Electrode molded Connector Shaft 2.4 mm Length 13.5 cm



EMI-1511
Silicone Monopolar 3M Cable



EMI-1574
Foot control pencil reusable with detachable 3M cable for 2.4mm Monopolar electrode



EMI-1573
Foot control pencil reusable with detachable 3M cable for 4 mm Monopolar electrode



GYNECOLOGY INSTRUMENTS



GRAVES VAGINAL SPECULUM

ITEM CODE	SIZE	LENGTH	WIDTH
EMI-3300	LARGE	150MM	35MM
	MEDIUM	95MM	35MM
	SMALL	35MM	20MM



PEDERSON VAGINAL SPECULUM

ITEM CODE	SIZE	LENGTH	WIDTH
EMI-3390	LARGE	150MM	25MM
	MEDIUM	100MM	25MM
	SMALL	70MM	14MM



GRAVES VAGINAL SPECULUM FOR LASER SURGERY

ITEM CODE	SIZE	LENGTH	WIDTH
EMI-3300-B	LARGE	150MM	35MM
	MEDIUM	95MM	35MM
	SMALL	35MM	20MM



COLLIN VAGINAL SPECULUM

ITEM CODE	SIZE	LENGTH	WIDTH
EMI-3460	LARGE	110MM	40MM
	MEDIUM	100MM	35MM
	SMALL	95MM	30MM



CUSCO VAGINAL SPECULUM

ITEM CODE	SIZE	LENGTH	WIDTH
EMI-3491	LARGE	95MM	35MM
	MEDIUM	85MM	35MM
	SMALL	75MM	32MM



GRAVES VAGINAL SPECULUM

ITEM CODE	SIZE	LENGTH	WIDTH
EMI-3300-S	LARGE	150MM	35MM
	MEDIUM	95MM	35MM
	SMALL	35MM	20MM



SILICONE BIPOLAR CABLES




EMI-1451



EMI-1450



EMI-1452



EMI-1471



EMI-1470



EMI-1490

DESCRIPTION

Silicone bipolar cable 3M, Instrument connector US 2-Pin plug generator connector valleylab, conmed, codman, and most bipolar generator.

EMI-1451

Silicone bipolar cable 3M, Instrument connector US 2-Pin plug generator connector 4mm banana Plug valleylab, conmed, codman, and most bipolar generator.

EMI-1450

Silicone bipolar cable 3M, Instrument connector US 2-Pin plug generator connector erbe, compatibility erbe siemens.

EMI-1452

Silicone bipolar cable 3M, Instrument connector european flat plug generator connector valleylab, conmed, codman, and most bipolar generator.

EMI-1471

Silicone bipolar cable 3M, Instrument connector european flat plug generator connector 4mm banana Plug valleylab, conmed, codman, and most bipolar generator.

EMI-1470

Silicone bipolar cable 3M, Instrument connector european flat plug generator connector erbe, compatibility erbe siemens.

EMI-1490

REUSABLE SCISSORS



EMI-5001-A-B

Dressing Scissors Length 14.5 cm Blunt / Blunt



EMI-5002-AB

Dressing Scissors Length 14.5 cm Sharp / Sharp Straight



EMI-5003

Dressing Scissors Length 18 cm Sharp / Sharp Straight



EMI-5004

Dressing Scissors Length 18 cm Blunt / Sharp Straight



EMI-5008

Metzenbaum Scissors Length 20 cm Curved



EMI-5009

Goldman Fox Scissors Curved Length 13 cm



EMI-5010

Goldman Fox Scissors Straight Length 13 cm



EMI-5011

Bush Umbilical Scissors Length 16 cm



EMI-5013

Iris Real Screw Scissors Length 11.5 cm Curved



EMI-5014

Iris Real Screw Scissors Length 11.5 cm Straight



EMI-5017

Iris Real Screw Scissors Length 12 cm Curved



EMI-5018

Iris Real Screw Scissors Length 12 cm Straight



EMI-5019

Spencer Stitch Scissor Length 16 cm



EMI-5020

Mayo Scissors Straight Length 17 cm



EMI-5021

Mayo Scissors Curved Length 17 cm



EMI-5027

Kelly Straight Scissor Length 16 cm



EMI-5028

Kelly Curved Scissor Length 16 cm



EMI-5029

Sims Utraine Scissors Curved Length 20 cm

REUSABLE FORCEPS



EMI-5045

Kelly Rankin Straight 16cm



EMI-5046

Kelly Rankin Curved 16cm




EMI-5047

Allis tissue Straight (3x4 teeth) 15 cm



EMI-5048

Luer Vulsellum Straight (1x1 teeth) 25 cm



EMI-5050

Sponge Holding Curved 25 cm



EMI-5051

Backhaus Towel Clip Forcep 9 CM



EMI-5053

Mayo Hegar Forcep 14, 16, 18 CM



EMI-5058

Halstead Mosquito Fine Forcep Curved 13 CM



EMI-5060

Spencer Wells Forcep Straight 13, 15, 18, 20 CM



EMI-5061

Ochsner Kocher Forcep Curved 14, 16, 18 CM



EMI-5062

Rochester Pean Forcep Curved 16 CM



EMI-5063

Kocher Forcep Straight 14, 16, 18 CM



EMI-5044

TC Iris Scissors Straight 12CM



EMI-5042

TC Dressing Scissors Curved, Sharp / Blunt, 14CM



EMI-5038

TC Metzenbaum Scissors Curved 18CM



EMI-5037

Halsey TC Forceps 13CM



EMI-5031

De baakey Forceps 16CM



EMI-5035

Kilner TC Forceps 18CM





ELECTROMED INSTRUMENTS (PVT) LTD.

ELECTROMED INSTRUMENTS (PVT) LTD.

HEAD OFFICE PAKISTAN:

-  +92 300 616 6061
+92 308 413 3335
-  +92 52 355 60 40
-  info@electromedsurgical.com
sales@electromedsurgical.com
-  +92 52 355 60 40
-  www.electromedsurgical.com
-  Mubarik Pura Road, Sialkot
(51310), Pakistan.

OFFICE AUSTRALIA:

-  +61421722813
-  2 Danehill place, Balga, Perth,
(6061) Western Australia.



CONTACT
QUICK SAVE

